

# Advancing Process Analytical Technologies

Exploring NIR for Robust High-Concentration  
Protein Analysis


**IFPAC Annual Meeting 2026**

Marija Iloska, Ph.D.

# Intros


## About Nirrin Technologies

- Founded in 2019
- HQ just north of Boston, MA



## Technical expertise


Tunable Near Infrared Lasers



Engineering and manufacturing




Analytics



# Talk agenda

1. The Shift to High-Concentration Biologics
2. UV-Vis Operating Regime at High Concentration
3. Spectroscopic Shift: Combination-Band NIR
4. From Univariate to Spectrally-Resolved Beer-Lambert
5. Real-World Deployment & UF/DF
6. Operational Impact & Future Directions



# High-Concentration Biologics Are Redefining Analytical Requirements

~30% of approved mAbs are formulated at  **$\geq 100$  mg/mL**


Many exceed **150–200 mg/mL** for subcutaneous delivery

Shawn Wang, Yifei Yan, Kin Ho, US FDA-approved therapeutic antibodies with high-concentration formulation: summaries and perspectives, *Antibody Therapeutics*,.

## Process & Analytical Implications

- Reduced measurement tolerance ( $\pm 1-2\%$  matters)
- Increased scattering & viscosity effects
- Greater UF/DF and formulation sensitivity
- Narrow operating window for UV-Vis at high concentration

# The UV-Vis Operating Regime at High Concentration



$$A(\lambda_{280}) = \epsilon_{protein}(\lambda_{280}) \cdot c_{protein} \cdot l$$

**Dynamic pathlength extends linear range**

**Classic conditions (25-50 mg/mL)**

- ~70um - 400um required pathlength

**New regime (150-250 mg/mL)**

- Required pathlength < 100um
- Measurement becomes geometry-, fouling- and scattering- sensitive

Figure: Iso-pathlength lines for typical IgG ( $\epsilon \cong 1.4 \text{ mg/mL}^{-1} \cdot \text{cm}^{-1}$ )

# The Spectroscopic Shift: From Electronic to Vibrational Absorption




Figure: UV iso-pathlength lines (black) for typical IgG ( $\epsilon_{280} \approx 1.4 \text{ mg/mL}^{-1} \cdot \text{cm}^{-1}$ ) and NIR (red  $\epsilon_{2290} \approx 3.3 \times 10^{-3} \text{ mg/mL}^{-1} \cdot \text{cm}^{-1}$ ) demonstrating Lower absorptivity preserves linear regime at high protein concentration.

UV (280 nm)	NIR (1200–2400 nm)
$\pi \rightarrow \pi^*$ electronic transitions	O–H, N–H, C–H vibrational modes
High extinction coefficient	Much lower absorptivity (400X lower $\epsilon$ )
Short pathlength required	mm scale fixed-pathlength viable
Aromatic residues only	Whole-molecule vibrational fingerprint

# Challenges of Conventional NIR in Aqueous Bioprocessing




Figure: Representative protein absorbance magnitudes for illustration

## UV-Vis: Analyte-dominated regime

- Protein absorption  $\gg$  water
- Direct Beer-Lambert quantitation
- High signal amplitude

## NIR: Background-dominated regime

- Water absorption dominates spectrum
- Protein appears as small perturbation
- Quantitation relies on differential sensitivity



# Enhancing Sensitivity in the NIR Region

## Overtone region (~1550-1800nm)


- Lower effective absorptivity
- Smaller concentration-dependent response

## Combination region (~2100-2350nm)

- Higher effective absorptivity
- Greater concentration-dependent response



# Pathlength Optimization in a Background-Dominated Regime



# Quantitative Protein Measurement in a Coupled Aqueous Matrix

## Physical Regime Established

- Linear absorbance regime at high concentration
- Combination-band region selected for maximal sensitivity
- Optimal pathlength identified for SNR

## Remaining Technical Challenges

- Overlapping spectral components (protein + water + buffer)
- Temperature-induced spectral variability
- Wavelength stability requirements
- Robust concentration estimation without extinction assumptions

# From Univariate to Spectrally Resolved Beer–Lambert

## Classical UV

$$A(\lambda_{280}) = \varepsilon_{protein}(\lambda_{280}) \cdot c_{protein} \cdot l$$

- Uses discrete or limited wavelengths for concentration inference
- Requires known absorptivity per wavelength
- Often assumes constant concentration during geometry adjustment
- Minimal spectral overlap

## Combination-band NIR

$$A(\lambda) = A_{measured}(\lambda) - A_{background}(\lambda)$$

$$A(\lambda) = \underbrace{\varepsilon_p(\lambda) l c_p}_{\text{protein}} + \underbrace{\varepsilon_b(\lambda) l c_b}_{\text{buffer}} + \underbrace{\sum_j \phi_j(\lambda) \beta_j}_{\text{nuisance base (path, temp, scatter)}}$$

**Full-spectrum fitting**

## Physically constrained multi-wavelength Beer-Lambert fitting

$$\mathbf{A}(\lambda_1 \dots \lambda_n) = \mathbf{E} \mathbf{c}$$

**Where:**

**A** = absorbance vector across wavelengths

**E** = matrix of basis spectra (water, protein, buffer)

**c** = concentration vector

$$\hat{\mathbf{c}} = \operatorname{argmin}_{\mathbf{c}} \|\mathbf{A}_{signal} - \mathbf{E} \mathbf{c}\|$$


# High-Precision Tunable Laser Spectroscopy (HPTLS)

## Architecture for Stable Spectral Basis Separation


HPTLS design element	Impact
Targeted combination-band region	Maximizes concentration-dependent contrast
Narrow-line tunable laser (sub-nm)	Resolves overlapping basis functions
Absolute wavelength referencing	Prevents basis misalignment
~20 $\mu$ AU repeatability	Preserves sub-milli-absorbance perturbations
Millimeter fixed pathlength	Maintains optimal SNR regime

**HPTLS enables physically constrained basis separation in a background-dominated regime**

# Transferable NIR Quantitation Across Bioprocess Formats




# Broad At-Line Dynamic Range IgG Quantitation Across Aqueous Environments



Reference concentration determined by orthogonal methods  
(SoloVPE, gravimetric, HPLC)

# Broad At-Line Dynamic Range IgG Quantitation Across Aqueous Environments




Reference concentration determined by orthogonal methods (SoloVPE, gravimetric, HPLC)

## Major takeaways


- Broad linear dynamic range (10 - 250 mg/mL)
- Fixed 1 mm pathlength (no dynamic pathlength adjustment)
- General IgG basis (NIST reference)
- Robust in varied aqueous formulations
- No dilution required

# Spectral Fits and Residual Analysis of 50 mg/mL Formulation


## Model basis contributions




## Measured vs Model Overlay



## Model Residual



# Integration of HPTLS Flow Cell in UFDF Skids




- Integrated in-line upstream of filter (pre-filtration position)
- Compatible with CIP/SIP, autoclave, or gamma-sterile tubing sets
- Modular flow-cell design for sterile bioprocess environments
- Fixed ~1 mm optical pathlength maintained during flow

# Real-Time Protein Quantitation During UFDF Concentration Ramps



## Real implementation

HPTLS flow cell connected to Akta Flux™ by Cytiva



## Key takeaways

- Demonstrated broad linear dynamic range (5 - 200 mg/mL) in flow
- General IgG basis (NIST reference)
- Robust in varied aqueous formulations

## Insights

- Fixed 1 mm pathlength maintains stable optical geometry during concentration ramps
- No pathlength adjustment required as absorbance increases

# Conclusions and Future Work

## What Combination-Band NIR Enables

- Fixed optical geometry across unit operations
- Broad linear range without dilution
- Transferable spectral basis across form factors
- Robust performance in background-dominated regimes
- Simplified PAT qualification pathway

## Next Steps Toward Scalable PAT

- Extend deployment to in-situ vessel monitoring
- Expanded therapeutic validation dataset
- Integration with control systems (closed-loop potential)
- Ongoing method validation & regulatory alignment